

Dandle ROO™ Neurodevelopmental Support System



Dandle ROO™

Dandle ROO™ was developed by clinicians to support the neurodevelopment of the preterm infant. It is an innovative nest system; complete with padded head boundary, gel pillow, shield and multipurpose roll.

Its unique design facilitates appropriate positioning and the soft, stretchable organic cotton surrounds and comforts the infant. Soft wings round the shoulders to encourage self-comforting behaviours whilst the lower wrap provides support and flexion to the lower extremities.

The multipurpose roll can fit in a pocket on the underside of the base to act as a prone roll or can be used in supine and side-lying positions to promote hip and pelvic flexion. A fitted gel pillow cushions and cradles the head and shoulders, helping to minimise head moulding while the detachable brim shields the infant's eyes from direct light.

In the Dandle ROO™ an infant can be easily moved with minimal distress for kangaroo care, weighing or examination.

Features

- Flexible aid supporting prone, supine and side line positioning
- Facilitates kangaroo care
- Encourages self-comforting behaviours
- Organic cotton, soft against the skin
- Multiuse it can be washed and reused
- Available in 3 sizes complete with corresponding gel pillow
- Optional laundry bag to help maintain integrity

Dandle ROO™ Neurodevelopmental Support System

TECHNICAL INFORMATION

Dandle ROO™ should be washed in accordance to the manufacturer's guidelines.

Made from organic cotton.



ORDERING INFORMATION

DESCRIPTION

Dandle ROO™

DESCRIPTION	PART NUMBER
Dandle ROO™ XS <1000g - complete (Case of 12)	DL 0951 B12
Dandle ROO™ XS <1000g - no gel pillow (Case of 12)	DL 0951 NG B12
Dandle ROO™ XS <1000g - packed in laundry bag (Case of 12)	DL 0951 LB B12
Dandle ROO™ S 1000-1800g - complete (Case of 12)	DL 0952 B12
Dandle ROO™ S 1000-1800g - no gel pillow (Case of 12)	DL 0952 NG B12
Dandle ROO™ S 1000-1800g - packed in laundry bag (Case of 12)	DL 0952 LB B12
Dandle ROO™ M 1800-2500g - complete (Case of 12)	DL 0953 B12
Dandle ROO™ M 1800-2500g - packed in laundry bag (Case of 12)	DL 0953 LB B12

Multi-Purpose Roll

Dandle ROO™ XS multi-purpose rolls (Case of 25)	DL 0951 ROLL B25
Dandle ROO™ S multi-purpose rolls (Case of 25)	DL 0952 ROLL B25
Dandle ROO™ M multi-purpose rolls (Case of 25)	DL 0953 ROLL B25

Gel Pillows

Dandle ROO™ XS gel pillow with cover 6" x 9" (Case of 12)	DL 0951 GP B12
Dandle ROO™ S gel pillow with cover 6" x 11" (Case of 12)	DL 0952 GP B12
Dandle ROO™ M gel pillow with cover 6" x 12" (Case of 12)	DL 0953 GP B12

Gel Pillow Covers

Dandle ROO™ XS gel pillow covers (Case of 12)	DL 0951 GPC B25
Dandle ROO™ S gel pillow covers (Case of 12)	DL 0952 GPC B25
Dandle ROO™ M gel pillow covers (Case of 12)	DL 0953 GPC B25

Laundry Bags

DandleLION™ laundry bags - 14" x 24" mesh (Case of 25)	DL 0981 B25
--	-------------

Above products are also available to purchase individually please contact Inspiration Healthcare directly for further details.

Nursing competencies are available to support the product range.

The details given in this leaflet are correct at time of going to press. The company reserves the right to improve the products shown.

Distributed by:



Connect Medizintechnik GmbH
Gspanngasse 4 | A-2130 Mistelbach
Tel: +43 664 168 68 08
Fax: +43 2572 20404
www.connect-medizintechnik.at



Inspiration Healthcare Ltd
Gildor House, West Street,
Earl Shilton, Leicester LE9 7EJ, UK

Tel: +44 (0) 1455 840 555

Fax: +44 (0) 1455 841 464

Email: info@inspiration-healthcare.com

Web: www.inspiration-healthcare.com

24 hour emergency hire and support call

+44 (0) 1455 840555

Dandle ROO™ is a trademark of Dandle LION Medical.