


Radiolucent NeoLead®

ECG Electrodes



N305



N306

Made in USA

FEATURES AND BENEFITS

- Pre-wired electrodes are nearly invisible on x-rays
- Hydrogel center
- NeoBond® Hydrocolloid border designed to:
 - prevent hydrogel from oozing or drying prematurely
 - provide sensor stability for better tracings
 - minimize slipping
- Multi-colored wires for easier identification
- N305 set of three leads
- N315 set of five leads for two channel, continuous monitoring
- Not made with natural rubber latex or plasticizer DEHP
- Sets of leads individually packaged

Micro NeoLead®

- Pre-wired electrodes are nearly invisible on x-rays
- Hydrogel center
- NeoBond® Hydrocolloid ends designed to:
 - provide sensor stability for better tracings
 - minimize slipping
 - reduce need for band or strap
- Multi-colored wires for easier identification
- Not made with natural rubber latex or plasticizer DEHP
- Sets of three individually packaged

Cat No.	Item	Size	Qty/Unit
N305	Radiolucent	1" Circle	20 sets/box
N306	Micro	1" x 3/8"	20 sets/box
N315	Pediatric	1" Circle	20 sets/box